

# Septal Button

Product Code: 890924

## Patient Information Leaflet (PIL)


This leaflet contains important information about your implant.


**You should speak to your doctor if you have any concerns or questions.**


All implantable medical devices have risks and benefits. Follow the advice from your healthcare team, even if it differs from the information contained in this leaflet. Please read this information carefully and keep it for future reference.

The name and number of your implant can be found on the patient implant card that has been provided by your doctor.

**OLYMPUS** Patient Implant Card

 \_\_\_\_\_

 \_\_\_\_\_

 \_\_\_\_\_

For more patient chart advice visit  
Patient information can be found at [www.olympus.com.au](http://www.olympus.com.au)

Device name \_\_\_\_\_

Model number \_\_\_\_\_

Batch / Lot number \_\_\_\_\_

Manufactured by

Olympus Medical Systems Corp. 2101 Mishawaka Dr.  
Irvine, CA 92614, USA. [www.olympus.com.au](http://www.olympus.com.au)

Olympus ACMI, Inc. 3800 Louisiana Avenue North Brooklyn Park  
MN 55445, USA. <http://medical.olympusamerica.com>

## Intended Use

The two cavities of your nose are separated by bone and cartilage. This is known as your septum, and it helps with airflow in the nasal passages.

**The Olympus Septal Button is used to close a perforation, or hole, that has developed in the septum.**

## Implant Description

The Septal Button is made from silicone. It can be trimmed to size by your doctor to best fit.

## Information for Safe Use

This implant can only be inserted by your doctor. The procedure may be performed in a clinic using local anaesthetic or can be combined with other surgical procedures.

You should follow your doctor's advice after your treatment. Discuss any questions or concerns with your doctor.

## Expected Lifetime and Follow Up

Your doctor will determine the length of time the Septal Button will be implanted based on your specific requirements.

## Adverse Effects and Reporting

Your doctor will provide information about the benefits, risks, and side effects of your treatment. Please talk to your doctor if you have concerns.

Possible adverse effects may include:

- Movement of the prosthesis
- Increased nosebleeds
- Crusting
- Pain within the nose
- Foreign body sensation
- Worsening of the hole in the septum

## Reporting adverse effects

If you wish to report any adverse effects you believe are a result of your treatment with this implant, please speak with your medical team or report the information to the Olympus Quality and Regulatory Department:

Email: [oaz-ra@olympus.com](mailto:oaz-ra@olympus.com)

Phone: 1300 132 992

You may also wish to contact the Therapeutic Goods Administration online at [www.tga.gov.au/reporting-problems](http://www.tga.gov.au/reporting-problems).

## Australian Sponsor

Olympus Australia Pty Ltd  
3 Acacia Place  
Notting Hill VIC 3168  
Australia

Phone: 1300 132 992

Website: <https://www.olympus.com.au/>

Patient Information:

<https://www.olympus.com.au/company/en/home/patient-information-leaflets.html>